Invacare | Platinum 9



Reduce cost with Invacare high flow concentrator

The Invacare Platinum 9 offers providers a value added and reliable global solution to treat patients with high flow requirements. Use as a stand alone device or couple together with HomeFill® Oxygen System for complete patient care.

While optimising oxygen concentration, Platinum 9 minimises costs associated to maintenance and delivery. It allows patients to stay at home comfortably and safely thanks to a series of integrated alarms.

Users will appreciate its simplicity of use combined with maximum security

Features and accessories







HomeFill® connection port

Technical data

Platinum 9

Platinum 9



365 mm



467 mm



670 mm



24 kg



230 V, 50 - 60 Hz





1 - 9 I/min



< 47 dB(A)



87-96%





Up to 1230m

ISO8359, MDD

SensO₂-System

O2-concentration



higher than 85% - green light -Normal



between 73% and 85% - yellow light -Check filter and flow rate



lower than 73% - red light & audible alarm -Call service technician

Service Interval

No provider routine preventative maintenance for 3 years, replace filters as required.

Ordering Information IRC9LXO2AWQ Platinum 9



